

# MAGNUM 2500

## MAGNETOTHERAPY | 52 PROGRAMS

PORTABLE MAGNETOTHERAPY, IDEAL FOR HOME TREATMENTS. MAGNUM 2500 IS RECOMMENDABLE FOR A BETTER RECOVERY FROM TRAUMAS AND FRACTURES, BUT ALSO TO GIVE RELIEF AND HELP ARTICULAR REGENERATION IN CHRONIC PATHOLOGIES SUCH AS ARTHROSIS AND ARTHRITIS.

### TECHNICAL SPECIFICATIONS

#### Display

Size of the visible area 2,6"

#### Channels

2 independents

#### Intensity can be emitted with flexible diffuser

- 160 Gauss on average for channel
- 320 Gauss on average on the whole (for 2 channels)

#### Frequency

5-200Hz

#### Duty Cycle

5-50%

#### Free programming

Possibility of adding up to 5 custom programs for each user

#### 1+1 Function

Possibility of running 2 different programs at the same time, using 1 channel each

#### Power supply

Mains

#### Size


100x160x35 mm

#### Weight

296 gr

### STANDARD EQUIPMENT

- 1 Bag
- 1 MAGNUM 2500 magnetotherapy device
- 1 Flexible solenoid 330x110 mm
- 1 Power supply unit
- 1 Operating manual



### PURCHASE OPTIONS

- G5438 MAGNUM 2500 with 1 flexible solenoid
- G5645 MAGNUM 2500 with 2 soft solenoids
- G5646 MAGNUM 2500 with 2 rigid solenoids



CE  
0476



LAST 10



FAVORITES



MULTI USER / 5 USERS



PROGRAMMABLE



RENTAL OPTION



1+1 MODE



RUN TIME / EDIT PARAMETERS



MAINS POWER SUPPLY

**GLOBUS**  
ITALIAN EXCELLENCE

# MAGNUM 2500

## MAGNETOTHERAPY | 52 PROGRAMS

### PROGRAM LIST

#### ARTHROSIS/ARTHRITIS

Knee arthrosis

Cervical arthrosis

Arthrotos

Arthrosis

Hand arthrosis

Hip arthrosis

Shoulder arthrosis

Column arthrosis

#### PAIN

Algodystrophy - Sudeck disease

Back pain

Joint pain

#### INFLAMMATION

Anti-inflammatory

Epicondylitis

Epitrochleytis

Carpal tunnel

Shoulder tendinitis

#### ORTHOPEDICS

Posto ACL surgery

Cartilage lesions

Knee prothesis

Hip prothesis

#### BONE PATHOLOGIES

Fractures

Scaphoid fracture

Wrist fracture

Shoulder Fracture

Collarbone fracture

Forearm fracture

Pelvic fracture

Tibial fracture

Femoral fracture

Delayed fracture consolidation

Fracture with plaster

Osteonecrosis

Osteoporosis

Bone edemas

#### TRAUMAS / EDEMAS / WOUNDS

Ankle sprain

Edemas- contusion

#### DEFAULT PROGRAMS WITH DIFFERENT FREQUENCY

Multifrequency program

Program 10Hz

Program 20Hz

Program 30Hz

Program 40Hz

Program 50Hz

Program 60Hz

Program 70Hz

Program 80Hz

Program 90Hz

Program 100Hz

Program 120Hz

Program 140Hz

Program 160Hz

Program 180Hz

Program 200Hz

