

MAGNUM 3500 PRO

PROFESSIONAL MAGNETOTHERAPY | 70 PROGRAMS

THE ONLY GLOBUS MAGNETOTHERAPY DEVICE WITH FOUR CHANNELS THAT ALLOWS TO TREAT DIFFERENT BODY PARTS OR SUBJECT AT THE SAME TIME. MAGNUM 3500 PRO IS EQUIPPED WITH A WIDE RANGE OF ACCESSORIES FOR ALL TREATMENT NEEDS.

TECHNICAL SPECIFICATIONS

Display

Size of the visible area 2.6"

Channels

4 independent channels

Intensity can be emitted with flexible diffuser

- 200 Gauss on average for channel
- 800 Gauss on average on the whole (for 4 channels)

Frequency

5-200Hz

Duty Cycle

5-50%

Free programming

Possibility of adding up to 10 custom programs for each user

2+2 function

Possibility of running 2 different programs at the same time, using 2 channels each

Power supply


Mains

Size | Weight

170x220x60 mm | 1152 gr

STANDARD EQUIPMENT

- 1 Suitcase
- 1 MAGNUM 3500 PRO magnetotherapy device
- 2 Flexible solenoids 330x110 mm
- 2 Soft solenoids 120x125 mm
- 1 Soft mat 4 250x250 mm
- 1 Power supply unit
- 1 Operating manual



CE
0476



- LAST 10
- FAVORITES
- MULTI USER / 20 USERS
- PROGRAMMABLE
- RENTAL OPTION
- 2+2 MODE
- RUN TIME / EDIT PARAMETERS
- MAINS POWER SUPPLY

GLOBUS
ITALIAN EXCELLENCE

MAGNUM 3500 PRO

PROFESSIONAL MAGNETOTHERAPY | 70 PROGRAMS

PROGRAM LIST

ARTHROSIS/ARTHRITIS

Knee arthrosis

Cervical arthrosis

Arthritis

Arthrosis

Hand arthrosis

Hip arthrosis

Shoulder arthrosis

Column arthrosis

PAIN

Algodystrophy

Back pain

Joint pain

INFLAMMATIONS

Anti-inflammatory

Tendinitis

Epicondylitis

Epitrochleytis

Carpal tunnel

Shoulder tendinitis

ORTHOPEDICS

Post ACL surgery

Cartilage lesions

Knee prosthesis

Hip prosthesis

Fusion

Osteotomy

BONE PATOLOGIES

Fractures

Scaphoid fracture

Wrist fracture

Shoulder fracture

Collarbone fracture

Forearm fracture

Pelvic fracture

Tibial fracture

Femoral fracture

Delayed fracture consolidation

Fracture with plaster

Osteobecrisis

Osteoporosis

Pseudoarthrosis

Bone edemas

TRAUMAS/EDEMAS/WOUNDS

Ankle sprain

Edemas - Contusion

Ulcers

NIGHT PROGRAMS

Night program 15Hz

Night program 50Hz

Night program 75Hz

DEFAULT PROGRAMS WITH DIFFERENT FREQUENCY

Multi-frequency program

Program 5Hz

Program 10Hz

Program 15Hz

Program 20Hz

Program 25Hz

Program 30Hz

Program 35Hz

Program 40Hz

Program 45Hz

Program 50Hz

Program 55Hz

Program 60Hz

Program 65Hz

Program 70Hz

Program 75Hz

Program 80Hz

Program 85Hz

Program 90Hz

Program 95Hz

Program 100Hz

Program 120Hz

Program 140Hz

Program 160Hz

Program 180Hz

Program 200Hz